




MOONS dummy (Order No. 08842)

Description

MOONS dummy to verify proper function of the MOONS setup.



In general, the MOONS dummy may be used only in combination with the cellasys MOONS setup, by qualified personnel from a research or healthcare institution.
Read the MOONS manual thoroughly and carefully follow the provided instructions and guidelines.

Values without MOONS dummy (open circuit):

	Value
Real part	NaN*
Imaginary part	NaN

MOONS dummy values (@ 10.000 Hz):

	Value
Real part	992 Ω +/- 20 Ω
Imaginary part	-121 Ω +/- 20 Ω

Technical data

Dimensions: 127,8 x 85,5 x 22,0 mm³

Weight: 48 g

Operating temperature: 0 °C to +80°C

R, C (parallel): 1 k Ω , 2,2nF

Calculated impedance: 981 Ω & -136 Ω @ 10 kHz

* NaN: Not a number

Liability / Copyright

All technical details are based on the state of the technology as of September 2026 and are subject to change without notice. No liability is assumed for pictures, descriptions, or any content of this document.

All descriptions, pictures, technical drawings, and all other illustrations are protected by copyright and, unless otherwise marked, are the property of cellasys know-how GmbH.

Any subsequent use requires prior written approval from cellasys know-how GmbH.

Contact

cellasys know-how GmbH
Illerstrasse 14
87758 Kronburg / Germany

www.cellasys.tech
info@cellasys.tech